

# Pure shift NMR: Applications



Mathias Nilsson

The University of Manchester

# Outline

---


The problem: spectral overlap in complex spectra

Pure shift NMR: achieving ultra-high resolution


Recovering coupling information when needed

Applications and assignment

# Pure shift NMR



Estradiol: 500 MHz  $^1\text{H}$  spectrum




Chemical shifts ( $\delta_{\text{H}}$ )

Homonuclear couplings ( $J_{\text{HH}}$ )

Overlap reduced by an order of magnitude


# 2D Pure shift NMR

Estradiol: 500 MHz TOCSY



Coupling information at pure shift resolution

# Pure shift NMR



Pure shift NMR suppresses multiplet structure to give a substantial improvement in resolution, comparable of going from a standard 500 MHz spectrometer to a hypothetical 5 GHz one

Adams, R. W. *EMagRes* 2007, 3 (4), 295–309.


Zangger, K. *Prog. Nucl. Magn. Reson. Spectrosc.* 2015, 86–87, 1–20.

Castañar, L. *Magn. Reson. Chem.* 2017, 55 (1), 47–53

M. Foroozandeh, G. A. Morris and M. Nilsson, *Chemistry-a European Journal*, 2018, **24**, 13988-14000.



J. A. Aguilar, S. Faulkner, M. Nilsson and G. A. Morris, *Angew. Chem. Int. Ed.*, 2010, **49**, 3901-3903.

# $^1\text{H}$ NMR of eye-drop formulation



Pure shift reveals minor components


# Pure shift with strong solvent signal



# Pure shift: revealing hidden peaks


- **Impurity** in valsartan formulation
- **Conventional HSQC/COSY**: impurity not evident (overlap/ similar spin systems)
- **Pure shift HSQC/COSY**: impurity becomes obvious.

Same acquisition time: COSY 2 min; HSQC 14 min



Pure shift reveals impurity correlations hidden by overlap

# Pure shift NMR – resolution vs coupling




✓ Resolution gained

× Coupling lost

Can we recover coupling information?

# Pure Shift NMR – 2DJ recovering lost coupling structure



2

3

4

5

6



7

8

## Phase sensitive 2DJ

- Full multiplets
- Moderate cost in sensitivity ( $\sqrt{2}$ )

# PSYCHEDELIC: individual J couplings





Selected-spin homonuclear J couplings – with pure shift resolution

# PSYCHEDELIC: multiple spins in one experiment


Colours = different selected spins

chenodeoxycholic acid



Multiple selected spins (not mutually coupled) in one experiment


# Pure shift TOCSY




2D correlations with pure shift resolution

# Pure shift MQF-COSY

Conventional




Pure shift



Mixture of 2R/2S isomers:  
hesperidin + naringin


2D pure shift enables isomer discrimination

# Pure shift real-time HSQC




Pure shift improves both resolution and sensitivity

# Pure shift HSQC in computer-assisted assignment




- 3-fucosyllactose (human milk oligosaccharide): **conventional** vs **pure shift** HSQC
- Automatic  $^1\text{H}/^{13}\text{C}$  assignment (CASPER): agrees with literature

$^1\text{H}$  assignment








$^{13}\text{C}$  assignment



One peak per  $\text{CH}_n$  correlation enables automated peak-picking


# Pure shift DOSY



PSYCHEiDOSY

Reduced overlap improves diffusion discrimination


# Analysis of complex mixtures



Sample + sequence




Peppermint oil



PSYCHE sequence

Pure shift provides ultra-high resolution  $^1\text{H}$  NMR spectra

# Selective 1D TOCSY



1D TOCSY resolution improved by pure shift

# Pure shift FESTA: resolving fluorinated spin systems by $^{19}\text{F}$ editing






Photo-degradation of betamethasone




FESTA selects the spin system; pure shift resolves every signal within it

# Decoupled pure shift NMR

(4S)-N-Boc-L-fluoroproline





(4R)-N-Boc-L-fluoroproline



$$J_{HH} = |0 - 25| \text{ Hz}$$


$$J_{HF} = |0 - 60| \text{ Hz}$$

<sup>1</sup>H NMR



**Pure shift**

Removes <sup>1</sup>H-<sup>1</sup>H couplings




**Fully decoupled**



Removes <sup>1</sup>H-<sup>1</sup>H couplings  
and <sup>1</sup>H-X couplings

Removing <sup>1</sup>H-X couplings reveals true pure shift resolution

# Broadband $^{19}\text{F}$ pure shift NMR




Mixture of toxic environmental fluorinated contaminants




Pure shift resolves chain-length defining  $\text{CF}_2$  groups in PFAS

# Broadband $^{31}\text{P}$ pure shift NMR



Isomers of  $(\text{PhNCO-SiMe}_3)_3\text{P}_7$



# Acknowledgements

---

## Manchester Group

Gareth A. Morris  
Laura Castañar  
Ralph W. Adams  
Adolfo Botana  
Mohammadali Foroozandeh  
Peter Kiraly  
Guilherme Dal Poggetto  
Pinelopi Moutzouri  
Maryam Kajeh  
Coral Mycroft  
Emma Gates  
Marshall Smith  
Hugo Rocha  
Nouran Hamed

## Collaborators

### *Industry*


Jonathan Bradley (J Matthey)  
Myron Johnson (J Matthey)  
Paul Bowyer (JEOL)  
Alexandra Shread (Bruker)

### *Academic*

Cláudio F. Tormena (Campinas)  
Anders Malmendal (Roskilde)  
Mikael Akke (Lund)  
Göran Widmalm (Stockholm)  
Davy Sinnaeve (Lille)  
Sebastian Meier (DTU)  
Jean-Nicolas Dumez (Nantes)  
Christina Thiele (Darmstadt)


*... and other collaborators  
(with thanks)*






# Resources from the Manchester NMR Methodology Group

---



**Pulse sequences**



**Software tools**



**Workshops & talks**



**Scan for resources**





[nmr.chemistry.manchester.ac.uk](http://nmr.chemistry.manchester.ac.uk)

# Pure shift methodology


## Pure shift key ideas




## Pure shift methods



# Basic pure shift sequence



Band selective




Sensitivity:

$\geq 100\%$

Broadband:

✗


Zangger-Sterk



0.5 – 10%

✓


PSYCHE



3 – 20%

✓

BIRD




$\approx 1\%$

✓

**ASR** enables pure shift; implementation determines sensitivity, bandwidth, and artefacts.


# Acquisition of pure shift – interferogram method

Acquire 'chunked' FID in increments



Interferogram

Real-time



Short chunks suppress J-evolution while preserving chemical shift evolution